

	<b>ANEXO INSTRUÇÃO DE TRABALHO</b>  <b>SISTEMA DE GESTÃO DA QUALIDADE</b>	Cód.: AIT-EXP-016
		Nº Rev.: 0
		Data Rev.: 24/05/2021
Assunto: TRANSPORT AND STORAGE POLICY Depto: EXPEDIÇÃO	Página: 1 / 1	

## TRANSPORT AND STORAGE POLICY

Wama Diagnóstica declares that temperature variations during the transport of products do not affect the quality of products, provided that the 96-hour interval is respected for products with refrigerator storage temperature (2 to 8°C) and 144-hour interval is respected for Imuno-Rápido line products (2 to 30°C).

In addition to real-time stability tests, product stability during transportation has been validated under non-recommended temperature stress conditions on the product packaging. The results of these studies show that the exposure of the products to temperatures different from the recommended ranges for a short period of time does not influence the performance of the kit until the expiration date of the kit.

This period of time is also included in the instruction of use and certificate of analysis of each product, determining the maximum shipping time outside the temperature described in the packaging of the kits.



Leise Daiana Martins  
Technical Manager  
CRBM-SP: 21,949